

Veterinary Screws, Drill Bits and Taps.

Quick Reference Chart













VET

PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES



Standard Screws, Drill Bits and Taps for Veterinary Use

SCREW DIAMETER (mm)										
Thread Diameter	1.5	1.5	2.0	2.4	2.7	3.5	4.5	5.5	4.0	6.5
Screw Type	Cortex									
						Cancellous				
	Full thread	Partial thread	Full thread	Full thread	Full thread	Full thread	Full thread	Full thread	Full thread	32 mm thread
Vet Part Number	VS102	VS105	VS202	VS203	VS205	VS302	VS402	–	VS403	VS603
Self-tapping	VS101	–	VS201	–	VS204	VS301	VS401	–	VS404	VS602
Non-self-tapping	–	–	–	–	–	–	–	VS501	–	–
DRILL BIT and TAP DIAMETER (mm)										
Drill Bit for Gliding Hole	1.5	1.5	2.0	2.4	2.7	3.5	4.5	5.5	–	4.5
GLIDING HOLE	⊗	⊗	⊗	⊗	⊗	⊗	⊗	⊗	⊗	⊗
THREADED HOLE										
Drill Bit for Threaded Hole	1.1	1.1	1.5	1.8	2.0	2.5	3.2	4.0	2.5	3.2
Tap	1.5	–	2.0	–	2.7	3.5	4.5	5.5	4.0	6.5
Catalog Number	311.15	–	311.19 311.21 (long)	–	311.26 311.28 (long)	311.32 311.33 (long)	311.46 311.48	311.55	311.34	311.66
Drive Type	1.5 mm Cruciform	T4 StarDrive™ Recess	2.0 mm Cruciform	2.4 mm Cruciform	2.5 mm Hex	2.5 mm Hex	3.5 mm Hex	3.5 mm Hex	2.5 mm Hex	3.5 mm Hex



Locking Screws and Drill Bits for Veterinary Use

SCREW DIAMETER (mm)							
Thread Diameter	1.5	2.0	2.4	2.7	3.5	4.0	5.0
Screw Type	Locking Head						
Vet Part Number Self-tapping	VS106	VS207	VS208	VS206	VS303	VS406	VS502
DRILL BIT and TAP DIAMETER (mm)							
Drill Bit for Gliding Hole	1.5	2.0	2.4	2.7	3.5	4.5	5.5
Drill Bit for Threaded Hole	1.1	1.5	1.8	2.0	2.8	3.2	4.3
Drive Type	T4 Stardrive Recess	T6 Stardrive Recess	T8 Stardrive Recess	T8 Stardrive Recess	T15 Stardrive Recess	T25 Stardrive Recess	T25 Stardrive Recess



DePuy Synthes
VET

PART OF THE FAMILY OF COMPANIES

DePuy Synthes Vet
1302 Wrights Lane East
West Chester, PA 19380
Telephone: (610) 719-5000
To order: (800) 523-0322

www.synthesvet.com

© DePuy Synthes Vet, a division of DOI 2018. All rights reserved.
088190-180228 03/18

